[Gebrauchsanleitung](#) | [Operating manual](#) | [Mode d'emploi](#) | [Instrucciones de manejo](#)

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Water jet pumps

Impressum

BRAND GMBH + CO KG
Otto-Schott-Str. 25
97877 Wertheim (Germany)

T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

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Contents

1	Terms of use	19
1.1	Hazard levels	19
1.2	Symbols used in the operating manual	19
1.3	Format	19
2	Safety	20
2.1	Operating Limitations.....	21
3	Functional elements.....	22
3.1	Application	22
4	Connection to water supply	23
5	Cleaning.....	25
6	Replace non-return valve.....	26
7	Technical data	27
8	Ordering Information	28
9	Warranty	30
10	Repairs.....	31
10.1	Sending for repair	31
11	Disposal	32

1 Terms of use

- Carefully read the operating manual before using the device for the first time.
- The operating manual is part of the device and must be kept in an easily accessible place.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find up-to-date versions of the operating manual on our website: www.brand.de.

1.1 Hazard levels

The following signal words identify possible hazards:

Signal word	Meaning
DANGER	Will lead to serious injury or death.
WARNING	May lead to serious injury or death.
CAUTION	May lead to minor or moderate injuries.
NOTICE	May lead to property damage.

1.2 Symbols used in the operating manual

Symbol	Meaning	Symbol	Meaning
	Danger area		Wear protective gloves
	Wear eye protection		

1.3 Format

Format	Meaning	Format	Meaning
1. Task	Indicates a task.	>	Indicates a condition.
a., b., c.	Indicates the individual steps of a task.	⇒	Indicates a result.

2 Safety

Please read carefully!

The laboratory device Water jet pump can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

1. Every user must read and understand this operating manual prior to using the device and observe these instructions during use.
2. Follow the general hazard instructions and safety regulations, e.g. wear protective clothing, eye protection and protective gloves. Observe special instructions for handling evacuated vessels, and if necessary use a protective shield.
3. Only use evacuated vessels which are definitely made for work under vacuum (e.g. desiccators and filter flasks).
4. Before evacuating, check all screw connections for absolute tightness.
5. In case of water pressure failure, dirty water must not be sucked back into the water mains. Therefore: Make sure that the water can always drain away freely. (Observe the DVGW* guidelines).
6. Always perform work in a manner that does not endanger users or other people.
7. Never use force! The fittings should be screwed by hand only (without using tools).
8. Use only original accessories and original replacement parts. Do not make any technical modifications. Do not dismantle the device any further than is described in the operating manual!

9. Always check that the device is in proper working condition before use. If there is a sign of a potential malfunction, immediately stop working with the device and follow cleaning and servicing instructions, see Cleaning, p. 25 Contact the manufacturer, if necessary.

*DVGW = German Technical and Scientific Association for Gas and Water

2.1 Operating Limitations

The user should always ascertain whether the device is suited for the particular application.


Some chemicals or their fumes may damage this device (swelling, embrittlement). Therefore: Do not use with chemicals which attack NBR, FKM, PP and PTFE.

Avoid exposure to temperatures exceeding 80 °C. Precision moulded parts of plastic material tend to distort at these temperatures.

The supply water pressure must be between 3 and 6 bar. At a lower pressure, the specified performance cannot be attained. A higher pressure may damage the pump.

The function of the non-return valve may be impaired through deposits or impurities, see Cleaning, p. 25 .

3 Functional elements


1	Hose nipple (water inlet connection)	2	Reducer
3	Sleeve nut 3/4"	4	Water inlet nozzle
5	Hose nipple (vacuum connection)		

3.1 Application

- Acceleration of filtering processes (vacuum filtration)
- Gentle distillation or sublimation (pressure reduction lowers the boiling or sublimation point)
- Drawing off liquids.
- Evacuation, e.g. of desiccators.


4 Connection to water supply

Threaded pipe connection

- Screw R $\frac{3}{4}$ " (2) sleeve nut tightly on water tap (1). A thread adapter for $\frac{1}{2}$ " with a gasket made out of rubber (NBR) is included. Reducers are available as accessories for $\frac{3}{8}$ " and M 22 x 1 (aerator thread), see Ordering Information, p. 28 .
- Screw the required thread adapter tightly on the water tap. Then connect the device. Do not use force.

Hose connection


- Unscrew the water inlet nozzle (with sleeve nut) from the water pump.
- Screw the water inlet nozzle (with hose nipple) tightly on the pump.


- Slide the water supply hose on the tube nipple till it fits tightly. Secure the water supply hose with a clip.

Performance check of the non-return valve

- a. Connect the filter pump to the water supply pipe (without vacuum line), open water inlet – the pump sucks.
- b. Close the water outlet by pressing your finger tip on it.
- c. No water should run out of the hose nipple now (vacuum arm).
- d. If there is, however, water leaking: renew the non-return valve, see Replace non-return valve, p. 26 .

Connection to vacuum line

- a. Before sliding on the vacuum tubing: Unscrew the hose nipple (vacuum arm) from instrument.
- b. Slide the vacuum tubing onto the loose hose nipple until it fits tightly.
- c. Re-fit the hose nipple on the filter pump.

NOTICE

Do not slide on the vacuum tubing as long as the filter pump is connected to the water tap.

5 Cleaning

Removal of internal sediments

- a. Unscrew the water inlet nozzle (thread or hose nipple).
- b. Unscrew the hose nipple for vacuum connection.

NOTICE

Do not dismantle the instrument any further!

- c. Immerse all parts into diluted hydrochloric acid (1:1) or a common deliming agent.
- d. Allow time for reaction, then rinse the parts thoroughly with water.
- e. Re-assemble all parts in the reverse order. Screw tight by hand. Do not use force.

6 Replace non-return valve

In case of water leakage during testing procedure of the non-return valve, the valve reed and PTFE sealing ring have to be replaced:


- a. Remove the hose nipple (vacuum arm) and O-ring.
- b. Pull out the PTFE sealing ring with a suitable tool (e.g. angled forceps).
- c. Remove the valve reed.
- d. Insert a new valve reed. Make sure that it is correctly placed in front of the valve opening.
- e. Slide in a new PTFE sealing ring right to the stop

New non-return-valve, see Ordering Information, p. 28 .


7 Technical data

Working fluid:	Water
Water supply line connection	R $\frac{3}{4}$ " ; with thread adapter R $\frac{1}{2}$ " ; With Olive: 10 - 12 mm external diameter
Connection to suction line	hose nipple, 6 - 9 mm external diam., with screw cap GL 14
Length of the filter pump:	approx. 210 mm (with $\frac{3}{4}$ " connection)
Weight	approx. 33 g (with $\frac{3}{4}$ " connection)
Material	Seal: NBR (nitrile butadiene rubber) O-ring: FKM (fluoro rubber material, replaceable) and NBR (permanently installed in the housing) Valve reed: FKM (fluoro rubber material) Sealing ring: PTFE All other parts: PP
Water consumption	approx. 220 litres/h at 3.5 bar (absolute)
Suction capacity	450 l/h (± 50 l/h) air against atmospheric pressure at supply pressure of 3.5 bar and a water temperature of 12 °C.
Ultimate pressure	at water temperature of 12 °C. ≤ 16 mbar
Water supply pressure	to reach the ultimate pressure: 3 - 6 bar (absolute)


Ultimate pressure as function of water temperature

Suction capacity as function of back pressure

Water consumption as function of supply pressure

8 Ordering Information

Description	Order No.
Water jet pump, complete with R $\frac{3}{4}$ " sleeve nut, reducer from R $\frac{3}{4}$ " to R $\frac{1}{2}$ ", hose connection (hose nipple 10 - 12 mm Ø) and vacuum connection (hose nipple 6 - 9 mm Ø) with valve plate made of FKM (fluoro rubber material)	159600

Accessories

Reducers

Description	Order No.
Reducer from R $\frac{3}{4}$ " to R $\frac{3}{8}$ " with seal, NBR	159665
Reducer from R $\frac{3}{4}$ " to M 22 x 1 (aerator thread), with seal, NBR	159670

Spares

Seal

Description	Order No.
Seal R $\frac{3}{4}$ "	159661
Seal R $\frac{3}{4}$ " (10 pcs.)	159664

Union nut

Description	Order No.
Sleeve nut R $\frac{3}{4}$ ", with Seal, NBR	159666

Reducer

Description	Order No.
Reducer from R $\frac{3}{4}$ " to R $\frac{1}{2}$ " with seal, NBR	159667

Hose nipple

Description	Order No.
Hose nipple for vacuum connection, 6 - 9 mm Ø, with sleeve nut and sealing, FKM	159669

Non-return valve

Description	Order No.
Non-return valve with valve reed out of FKM, complete	159668

9 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the instrument or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

USA and Canada:

Find more warranty information on www.brandtech.com.

10 Repairs

10.1 Sending for repair

NOTICE

Transporting hazardous materials without approval is prohibited by law.

Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

Outside USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

Contact addresses

Germany:

BRAND GMBH + CO KG
Otto-Schott-Straße 25
97877 Wertheim (Germany)
T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T +91 22 42957790
F +91 22 42957791
info@brand.co.in
www.brand.co.in

USA and Canada:

BrandTech® Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)
T +1-860-767 2562
F +1 - 860 - 767 2563
info@brandtech.com
www.brandtech.com

China:

BRAND (Shanghai) Trading Co., Ltd.
Rm 201-202, North Tower,
No. 199 Kaibin Rd, Xuhui District, Shanghai
Shanghai 200030 (P.R. China)
T +86 21 6422 2318
F +86 21 6422 2268
info@brand.com.cn
www.brand.cn.com

11 Disposal

Before disposal, observe the relevant national disposal regulations and ensure that the product is disposed of properly.

Subject to technical changes, errors, and misprints.

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