



Seal-R-film™

The stretchable universal sealing film

BRAND. For lab. For life.®

- + Protects against evaporation, leakage, and contamination
- + Easy sealing with 300% elasticity
- + Lot identification of each roll

Sealing film Seal-R-film™

No more spilling, evaporation, and contamination!

Seal-R-film™ sealing film seals beakers, flasks, and other vessels. The sealing film can be stretched to three times its original length. This means that Seal-R-film™ always fits, even with irregular shapes. For traceability, each roll has a lot identification.

With the BRAND cutter for sealing film you can cut rolls with widths of 50 and 100 mm efficiently and economically and store them cleanly.

Features:

- + Thickness: 0.14 +/- 0.02mm
- + Toxicity: non-toxic
- + Melting point: > 66°C (151°F)
- + Boiling point: > 288°C (550°F)
- + Continuous service temperature: -20°C to 46°C
- + Elasticity (25°C): 300%
- + Tear strength (25°C): 1.5 MPa

Gas permeability in 24h at room temperature and a pressure of 0.1 MPa:

CO₂ (carbon dioxide) < 1000 cm³/m²

O₂ (oxygen) < 350 cm³/m²


Water vapor permeability < 0.8 g/m²

Easy to use

Simply take Seal-R-film™ from the roll, remove the paper, and pull it over the area to be sealed while stretching the film uniformly.



All product information
can be found at
shop.brand.de




Ordering information

Seal-R-film™

Universal sealing film

Width x Length	Pack of roll	Cat. No.
100 mm x 38 m	12	701620
	1	701621
50 mm x 75 m	24	701622
	1	701623
100 mm x 75 m	12	701624
	1	701625
500 mm x 15 m	6	701626
	1	701627



Cutter for sealing film

Suitable for 50 mm and 100 mm wide rolls.

Description	Pack of	Cat. No.
Cutter	1	701650

Seal-R-film™ applications

- + Protecting sensitive samples from contamination
- + Sealing containers such as flasks, beakers, cuvettes, petri dishes and irregularly shaped containers
- + Covering plugs to ensure a very secure seal
- + Protecting from humidity for long-term storage

Seal-R-film™ chemical resistance

The sealing film is resistant up to 48 hours against many polar substances, e.g., saline solutions, inorganic acids and alkaline solutions. After this period embrittlement may occur.

Effects of 48-hour exposure at 30°C, relative humidity 80%.

Medium	
Acids	
Hydrochloric acid (HCl) 36,5%	resistant
Sulfuric acid (H ₂ SO ₄) 98%	resistant
Nitric acid (HNO ₃) 95%	resistant*
Hydrogen peroxide (H ₂ O ₂) 30%	resistant
Alkalis	
Sodium hydroxide (NaOH) 22%	resistant
Ammonia (NH ₃) 28%	resistant
Saline solutions	
Sodium chloride (NaCl) 20%	resistant
Potassium permanganate (KMnO ₄) 5%	resistant*
Iodine solution 0.1 mol/l	resistant*
Organic solvents	
Methanol	resistant
Ethanol	resistant
Isopropanol	resistant

* Brown coloration

BRAND GMBH + CO KG

P.O. Box 1155 | 97861 Wertheim | Germany

T +49 9342 808 0 | F +49 9342 808 98000 | info@brand.de | www.brand.de



BRAND. For lab. For life.®

BRAND®, BRAND. For lab. For life.®, as well as the BRAND figurative mark are registered trademarks or trademarks of BRAND GMBH + CO KG, Germany. All other trademarks mentioned or depicted here are the property of the respective owners.

Our technical literature is intended to inform and advise our customers. However, the validity of general empirical values, and of results obtained under test conditions, for specific applications depends on many factors beyond our control. Please appreciate, therefore, that no claims can be derived from our advice. The user is responsible for checking the appropriateness of the product for any particular application.

California Residents: For more information concerning California Proposition 65, please refer to www.brand.de/calprop65.

Subject to technical modification without notice. Errors excepted.

994140 © 2023 BRAND GMBH + CO KG | Printed in Germany | 0423



Find accessories and replacement parts, user manuals, test instructions (SOP) and product videos at shop.brand.de



Further information on products and applications can be found on our YouTube channel: [mylabBRAND](https://www.youtube.com/mylabBRAND)



BRAND (Shanghai) Trading Co., Ltd.
Shanghai, China

T +86 21 6422 2318
info@brand.com.cn
china.brand.com.cn

BRAND Scientific Equipment Pvt. Ltd.
Mumbai, India

T +91 22 42957790
customersupport@brand.co.in
www.brand.co.in

BrandTech® Scientific, Inc.
Essex, CT. United States of America

T +1 860 767 2562
info@brandtech.com
www.brandtech.com