



TipStack

Less waste – space-saving design.
Stacking towers are back!

BRAND. For lab. For life.®

TipStack

A tip tower containing 5 filled tip-trays and a TipBox constitute the space-saving refill system for 20 µl, 200 µl and 1000 µl tips. Tightly sealing spacers prevent the tips from getting stuck together, and ensure them to be free from DNA, RNases, endotoxins and ATP. The sterile TipStacks (BIO-CERT® quality) are supplied with a transfer aid for contamination-free use in a previously sterilized TipBox.

Your benefits:

- + All components are recyclable
- + Reduced amount of waste
- + Sterilizable and reusable TipBox
- + Space-saving design
- + With transfer aid for sterile tips
- + BIO-CERT® LIQUID HANDLING QUALITY




The tip system

The space-saving refill system, consisting of a tip tower with five filled tip-trays and a TipBox, offers an efficient solution for storing and refilling pipette tips in sizes 20 µl, 200 µl, and 1000 µl. All components of the system are recyclable, thus contributing to a reduction in waste. The TipBox itself is sterilizable and reusable, reducing both costs and environmental impact. Ideal for laboratories that value efficiency, sustainability, and the highest quality.



TipStack



Transfer aid



TipBox



All product information can be found at shop.brand.de

Sterile handling

Ordering information

TipStack, pipette tips, racked, BIO-CERT® LIQUID HANDLING QUALITY

	Volume	Pack of	non-steril Cat. No.	sterile Cat. No.
	0.5 – 20 µl	960 pipette tips (2 x 5 TipRacks, 96 each)	732244	732264
	2 – 200 µl	960 pipette tips (2 x 5 TipRacks, 96 each)	732248	732268
	50 – 1000 µl	960 pipette tips (2 x 5 TipRacks, 96 each)	732252	732272



EcoVadis silver rating 2025:
BRAND is among the top 15%
of comparable companies.

BRAND®, BRAND. For lab. For life.®, as well as the figurative marks depicted here and the BRAND figurative mark are registered trademarks or trademarks of BRAND GMBH + CO KG, Germany. All other trademarks mentioned or depicted here are the property of the respective owners. Subject to technical modification without notice. Errors excepted.

BRAND GMBH + CO KG | P.O. Box 1155 | 97861 Wertheim | Germany
T +49 9342 808-0 | info@brand.de | www.brand.de

