

The Transferpette® S-8 and S-12 multi-channel pipettes

Work in microplates and PCR plates is efficient with the multi-channel models of the Transferpette® S. These 8- and 12-channel pipettes cover the volume range

from 0.5 µl to 300 µl. They provide all the advantages of the single-channel models, plus additional features that make work and maintenance easier.



Adapts to you

Pipetting unit can be rotated 360° for individualized comfort

Low forces

Easy mounting and ejection of tips due to O-rings

Easy cleaning and maintenance

Individual nose cones and seals can be easily removed for cleaning or replacement



Narrow centrifuge tubes



Sample preparation



Working with PCR plates



Easy maintenance