

At a glance: Advantages of the Transferpette® S

Whether you are using the single-channel pipette or the 8- and 12-channel models, the Transferpette® S is always the first choice for simple and efficient pipetting in demanding applications.

The single-channel pipettes cover a wide volume range from 0.1 µl to 10 ml. Pipettes with small nominal volumes provide high precision for small-volume applications such as those found in molecular biology, e.g. for pipetting enzymes.

Large nominal volumes are used for water and environmental analysis, for example, where durability is critical. For standardized applications with specific volumes, we offer economical fixed volume models.

Genuine one-hand operation

Set the volume, pipette and eject the tip – all without changing your grip

Volume-change protection

Secures the set volume from accidental adjustment

Ergonomic

Low operating forces and short stroke

Easy Calibration technology

Simple adjustment, without tools

MyProduct

Scan the QR code and request certificates and documents for your individual product online

Transferpette® S,
fixed volume

Transferpette® S,
adjustable volume

Slim shaft

Easy access to narrow vessels



Easy Calibration technology



Volume-change protection



Integrated shaft coupling



Product documentation