


# 1 – 50 µl Pipette tips

## 732006 | 732026




Dimensions in mm

<b>Volume</b>	1 – 50 µl
<b>Material</b>	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
<b>Color</b>	transparent
<b>Packaging</b>	bulk
<b>Quality</b>	Produced under cleanroom conditions acc. to ISO 14 644-1 class 8



# 1 – 50 µl Pipette tips

## 732106 | 732206



Dimensions in mm


<b>Volume</b>	1 – 50 µl
<b>Material</b>	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
<b>Color</b>	transparent
<b>Packaging</b>	racked
<b>Quality</b>	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **



\* human and bacterial DNA, \*\* according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

# 1 – 50 µl Pipette tips

## 732126 | 732226



Dimensions in mm

<b>Volume</b>	1 – 50 µl
<b>Material</b>	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
<b>Color</b>	transparent
<b>Packaging</b>	racked
<b>Quality</b>	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of ATP



\* human and bacterial DNA, \*\* according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml