


50 – 1000 µl Pipette tips 732012 | 732032




Dimensions in mm

Volume	50 – 1000 µl
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	blue
Packaging	bulk
Quality	Produced under cleanroom conditions



50 – 1000 µl Pipette tips

732112 | 732212



Dimensions in mm


Volume	50 – 1000 µl
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

50 – 1000 µl Pipette tips

732132 | 732232



Dimensions in mm

Volume	50 – 1000 µl
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of ATP



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml