


# UV-Cuvette, micro

## 759200, 759210, 759215

PRODUCT  
INFORMATION



Dimensions in mm


<b>Filling volume min.</b>	70 $\mu$ l
<b>Filling volume max.</b>	850 $\mu$ l
<b>Light path</b>	10 mm
<b>Center height</b>	8.5 mm
<b>Material</b>	UV-polymer

Subject to technical modification without notice. Errors excepted.

# UV-Cuvette, micro

## 759220, 759230, 759235

PRODUCT  
INFORMATION



Dimensions in mm

<b>Filling volume min.</b>	70 $\mu$ l
<b>Filling volume max.</b>	550 $\mu$ l
<b>Light path</b>	10 mm
<b>Center hight</b>	15 mm
<b>Material</b>	UV-polymer

Subject to technical modification without notice. Errors excepted.