

Bulb pipettes

Capacity [ml]	Color Code (ISO 1769)	Pattern	Length mm [± 10 mm]	Error limit [± ml]		
				BLAUBRAND®	BLAUBRAND® USP	SILBERBRAND
0.5	2x black	straight	300	0.005	0.005	0.007
1	blue	straight	300	0.008	0.006	0.010
2	orange	bulb	330	0.010	0.006	0.015
2.5	–	bulb	330	0.010	–	–
3	black	bulb	330	0.010	0.010	–
4	2x red	bulb	400	0.015	0.010	–
5	white	bulb	400	0.015	0.010	0.025
6	2x orange	bulb	400	0.015	0.010	–
7	2x green	bulb	400	0.015	0.010	–
8	blue	bulb	440	0.02	0.02	–
9	black	bulb	440	0.02	0.02	–
10	red	bulb	440	0.02	0.02	0.03
15	green	bulb	510	0.03	0.03	–
20	yellow	bulb	510	0.03	0.03	0.04
25	blue	bulb	520	0.03	0.03	0.04
30	black	bulb	520	0.03	0.03	–
40	white	bulb	540	0.05	0.05	–
50	red	bulb	540	0.05	0.05	0.07
100	yellow	bulb	585 (± 15 mm)	0.08	0.08	0.12

