

Cylinder

Graduated cylinder: tall form

Capacity [ml]	Subdivision [ml]	Height [mm]	Outer Ø [mm]	Width across corners x width across flats [mm]	Error limit [± ml]		
					BLAUBRAND®	BLAUBRAND® USP	SILBERBRAND
5	0.1	115	12	50 x 44	0.05	0.05	0.08
10	0.2	140	14	50 x 44	0.10	0.10	0.15
25	0.5	170	21	60 x 53	0.25	0.17	0.4
50	1	200	25	63 x 56	0.5	0.25	0.8
100	1	260	29	76 x 67	0.5	0.5	0.8
250	2	335	39	102 x 90	1.0	1.0	1.5
500	5	365	53	113 x 100	2.5	2.0	4
1000	10	465	65	141 x 125	5	3.0	8
2000	20	505	85	158 x 140	10	6.0	15

Graduated cylinder: low form

Capacity [ml]	Subdivision [ml]	Height [mm]	Outer Ø [mm]	Width across corners x width across flats [mm]	Error limit [± ml]	
					SILBERBRAND	
10	1	90	21	60 x 53	0.3	
25	1	115	25	63 x 56	0.5	
50	2	145	29	76 x 67	1	
100	2	165	39	102 x 90	1	
250	5	195	54	113 x 100	2	
500	10	250	65	141 x 125	5	
1000	20	285	85	158 x 140	10	
2000	50	340	105	175 x 156	20	

Mixing cylinder: tall form

Capacity [ml]	Subdivision [ml]	Height [mm]	Outer Ø [mm]	Width across corners x width across flats [mm]	Error limit [± ml]	
					BLAUBRAND®	SILBERBRAND
10	0.2	160	14	50 x 44	0.10	0.15
25	0.5	190	21	60 x 53	0.25	0.4
50	1	220	25	63 x 56	0.5	0.8
100	1	285	29	76 x 67	0.5	0.8
250	2	350	39	102 x 90	1.0	1.5
500	5	395	53	113 x 100	2.5	4
1000	10	500	65	141 x 125	5	8

