

Graduated pipettes

Capacity [ml]	Subdivision	Color Code (ISO 1769)	Length mm [± 10 mm]	Error limit [± ml]		
				BLAUBRAND®	BLAUBRAND® USP	SILBERBRAND
0.5	0.01	2x yellow	360	0.006	–	0.008
1	0.01	yellow	360	0.007	0.007	0.010
1	0.1	red	360	0.007	–	0.010
2	0.01	2x white	360	0.010	–	0.015
2	0.02	black	360	0.010	0.010	0.015
2	0.1	green	360	0.010	–	0.015
5	0.05	red	360	0.030	0.020	0.05
5	0.1	blue	360	0.030	–	0.05
10	0.1	orange	360	0.05	0.03	0.08
20	0.1	2x yellow	360	0.1	–	0.15
25	0.1	white	450	0.1	–	0.15
50	0.5	–	450	0.2	–	–

